


ALP TITANIUM PROXIMAL HUMERUS PLATE

Proh-Lock



PERFECT BONE FIT



ALP TITANIUM DISTAL LATERAL & MEDIAL HUMERUS PLATE


TELIOS

TRAUFIX
READY TO RUN

ALP TITANIUM DISTAL LATERAL & MEDIAL HUMERUS PLATE | TELIOS

ANATOMICAL STRUCTURE

OFFSET BOLTS ON THE BODY OF THE PLATE TO AVOID BOLT COLLISION



LOW DISTAL PROFILE **0.70mm**
(Lateral & Medial Plate)

8 HOLES OF LOCKING
VARIABLE ANGLE
(Distal lateral Telios)

Distal Lateral Telios

Left (Code 145) / Right (Code 144) Version

3, 5, 7 and 9 Shaft holes

Distal Medial Telios

Left (Code 147) / Right (Code 146) Version

2, 4, 6 and 8 Shaft holes

2.5mm TITANIUM CORTEX LOCKING SCREW, TORX DRIVE (Code 155X, Length 08mm-30mm)
3.5mm TITANIUM CORTEX LOCKING SCREW (Code 106, Length 12mm-60mm)
3.5mm TITANIUM CANCELLOUS LOCKING SCREW (Code 107, Length 12mm-60mm)
3.5mm TITANIUM CORTEX SCREW (Code 112, Length 12mm-60mm)

ALP TITANIUM PROXIMAL HUMERUS PLATE | PROH-LOCK

LOW PROXIMAL PROFILE **2.5mm**

100% ANATOMICAL

LOW DISTAL PROFILE **3.5mm** Proh-Lock



LOW DISTAL PROFILE **4mm** Proh-Lock Large



ROUNDED PROFILE PLATE

9 PROXIMAL POLYAXIAL HOLES OF
LOCKING VARIABLE ANGLE $\pm 15^\circ$

15° 15°

Proh-Lock

Left (Code 130L) / Right (Code 130R) Version

3 and 5 Shaft holes

Proh-Lock Large


Left (129L) / Right (Code 129R) Version

4, 6, 8 and 10 Shaft holes

3.5mm TITANIUM CORTEX LOCKING SCREW (Code 106, Length 12mm-60mm)
3.5mm TITANIUM CANCELLOUS LOCKING SCREW (Code 107, Length 12mm-60mm)
3.5mm TITANIUM CORTEX SCREW (Code 112, Length 12mm-60mm)


DISTAL LATERAL & MEDIAL TELIOS

Specifications



Total Lengths Chart

PLATE	A) LATERAL TELIOS	PLATE	B) MEDIAL TELIOS
HOLES	LT (mm)	HOLES	LT (mm)
3 Holes	107.00	2 Holes	106.00
5 Holes	135.00	4 Holes	134.00
7 Holes	163.00	6 Holes	162.00
9 Holes	191.00	8 Holes	190.00



Titanium Alloy 6Al4V ELI ASTM F136

Special for surgical implants

Design by **TRAUFIX**
READY TO RUN



MADE IN MEXICO



*Refer to 510(k) for cleared devices. Some pictures of devices in this banner are not 510k cleared.



510 (k) K201510
*Consult products

ALP MODULAR SYSTEM | ALL IN ONE KIT

PLATES AND INSTRUMENTS IN ONE KIT | EASY & COMFORTABLE

ALP SMALL FRAGMENTS SET


QTY. EQUIPMENT

- 1 Eccentric neutral Drill guide Ø 2.5mm
- 2 Ø 2.5 / 3.5 Double Drill guide
- 1 Hexagon screwdriver tip with AO anchor Ø2.5mm
- 1 Screw extractor tip (with cone tip) Ø2.5mm
- 1 Countersink tip Ø6mm
- 1 Reamer Tip (Punch) for Small Fragments Ø 6.5 mm
- 1 T-Handle with quick coupling (AO)
- 1 Tap point Ø3.5mm
- 1 Handle with torque wrench AO anchor Ø1.5 NM
- 1 Depth gauge stainless steel Ø60mm
- 2 Guide wire with threaded tip Ø1.1mm
- 2 Kirschner Wire Ø1.6mm
- 1 Screwdriver with hexagon tip for 3.5/4.0 w/holder screw Ø2.5mm
- 1 Drill bit with AO anchor Ø2.5 x 150mm
- 1 Drill bit with AO anchor Ø2.7 x 150mm
- 1 Drill bit with AO anchor Ø3.5 x 150mm
- 2 Drill guide / sleeve Ø2.7mm

ALP SUPPORT EQUIPMENT FOR SMALL FRAGMENTS

QTY. EQUIPMENT

- 2 Hohmann spacers (pair) 8/22
- 2 Mini lane
- 2 Multitoothed tweezers
- 2 Field type clamp 13.5mm
- 1 Kirschner cutter f/pin kirschner
- 1 Volkman spoon 17cm
- 1 Gouge forceps 20 cm
- 1 Triscators for small fragments
- 1 Osteotomo mini lambotte 8mm
- 1 Osteotomo mini lambotte 12mm
- 1 Periosteal elevator 9 mm
- 2 Verbrugge clamp 17.5cm
- 1 Hammer 200 G
- 1 Field type clamp 20cm
- 1 16 cm multitoothed tweezers
- 1 Stop end clamp 120 cm
- 2 Hexagon screwdriver tip with AO T8 anchor Ø2.5mm
- 1 Handle with torque AO green 0.8 NM
- 1 Osteotome mini-lamb Ø3 x 125mm
- 1 Osteotome mini-lamb Ø4 x 125mm
- 1 Osteotome mini-lamb Ø6 x 125mm
- 1 Osteotome mini-lamb Ø8 x 125mm
- 1 Watchmaker handle with blue AO input Ø2.5mm
- 2 Drill guide / sleeve Ø2.0 / 2.5mm
- 2 Drill bit with AO anchor Ø2.0 x 125mm
- 1 Lionter gauge 50mm
- 1 Ø 2.0/2.7mm Double Drill guide
- 1 Adson brown forceps
- 2 Mini-Hohmann spacer Ø 6mm
- 2 Kirschner Wire Ø1.1 X 100mm



internationalsales@fixier.com.mx
Mobile: (+52) 55 8106 3240

www.traufix.com

 traufix.mx

 traufix