

FIBULE

ALP TITANIUM DISTAL FIBULA PLATE FIBULE



PERFECT BONE FIT


TIDIS II

ALP TITANIUM DISTAL PLATE TIDIS II



TRAUFIX
READY TO RUN

ALP TITANIUM DISTAL TIBIA PLATE TIDIS II



ROUNDED PROFILE PLATE
TO REDUCE THE RISK OF SOFT TISSUE IRRITATION


LOW DISTAL PROFILE **1.6mm**



SEVEN POLYAXIAL HOLES ($\pm 15^\circ$) AND TWO MONOAXIALS
TO AVOID JOINT DAMAGE TO THE FIBULA

DISTAL HOLE TOWARDS THE TIBIAL MALLEOLUS WITH
30° DEGREE INCLINATION + 15° DEGREE POLYAXIAL

SUPERFICIAL LOW CONTACT
NOTCHES AT THE TOP OF THE PLATE
ALLOW BLOW VESSEL CIRCULATION



Distal Tibia Plate Tidis II

Left(Code 165) / Right(Code 164) Version

5, 7, 9, 11 and 13 Shaft holes

3.5mm TITANIUM CORTEX LOCKING SCREW (Code 106, Length 12mm-60mm)
3.5mm TITANIUM CANCELLOUS LOCKING SCREW (Code 107, Length 12mm-60mm)
3.5mm TITANIUM CORTEX SCREW (Code 112, Length 12mm-60mm)

ALP TITANIUM DISTAL FIBULA PLATE FIBULE



FIVE DISTAL POLYAXIAL HOLES THAT ACCEPT 2.5mm AND 2.7mm
LOCKING SCREWS

COMBINED HOLES

ROUNDED EDGES TO MINIMISE SOFT ISSUE IRRITATION

ANATOMICAL STRUCTURE


Distal Fibula Plate Fibule

Left(Code 234) / Right(Code 233) Version

3, 4, 5, 6, 7, 9 and 11 Shaft holes

2.5mm TITANIUM CORTEX LOCKING SCREW (Code 155, Length 08mm-30mm)
2.7mm TITANIUM CORTEX LOCKING SCREW (Code 170, Length 06mm-60mm)
3.5mm TITANIUM CORTEX LOCKING SCREW (Code 106, Length 12mm-60mm)
3.5mm TITANIUM CORTEX SCREW (Code 112, Length 12mm-60mm)


Plate design is anatomical and low profile, both
distally and diaphyseal.




DISTAL TIBIA PLATE TIDIS II

Total Lengths Chart

PLATE	TIDIS II
HOLES	LT (mm)
5 Holes	121.00
7 Holes	147.00
9 Holes	173.00
11 Holes	199.00
13 Holes	225.00



Specifications




FIBULA PLATE FIBULE

Total Lengths Chart

PLATE	FIBULE
HOLES	LT (mm)
3 Holes	73.00
4 Holes	86.00
5 Holes	99.00
6 Holes	114.00
7 Holes	127.00
9 Holes	153.00
11 Holes	179.00

Specifications



Titanium Alloy 6Al4V ELI ASTM F136

Special for surgical implants

Design by **TRAUFIX**
READY TO RUN



MADE IN MEXICO



*Refer to 510(k) for cleared devices. Some pictures of devices in this banner are not 510k cleared.



510 (k) K201510
*Consult products

ALP MODULAR SYSTEM | ALL IN ONE KIT

PLATES AND INSTRUMENTS IN ONE KIT | EASY & COMFORTABLE

ALP SMALL FRAGMENTS SET


QTY. EQUIPMENT

- 1 Eccentric neutral Drill guide Ø 2.5mm
- 2 Ø 2.5 / 3.5 Double Drill guide
- 1 Hexagon screwdriver tip with AO anchor Ø2.5mm
- 1 Screw extractor tip (with cone tip) Ø2.5mm
- 1 Countersink tip Ø6mm
- 1 Reamer Tip (Punch) for Small Fragments Ø 6.5 mm
- 1 T-Handle with quick coupling (AO)
- 1 Tap point Ø3.5mm
- 1 Handle with torque wrench AO anchor Ø1.5 NM
- 1 Depth gauge stainless steel Ø60mm
- 2 Guide wire with threaded tip Ø1.1mm
- 2 Kirschner Wire Ø1.6mm
- 1 Screwdriver with hexagon tip for 3.5/4.0 w/holder screw Ø2.5mm
- 1 Drill bit with AO anchor Ø2.5 x 150mm
- 1 Drill bit with AO anchor Ø2.7 x 150mm
- 1 Drill bit with AO anchor Ø3.5 x 150mm
- 2 Drill guide / sleeve Ø2.7mm

ALP SUPPORT EQUIPMENT FOR SMALL FRAGMENTS

QTY. EQUIPMENT

- 2 Hohmann spacers (pair) 8/22
- 2 Mini lane
- 2 Multitoothed tweezers
- 2 Field type clamp 13.5mm
- 1 Kirschner cutter f/pin kirschner
- 1 Volkman spoon 17cm
- 1 Gouge forceps 20 cm
- 1 Triscators for small fragments
- 1 Osteotome mini lambotte 8mm
- 1 Osteotome mini lambotte 12mm
- 1 Periosteal elevator 9 mm
- 2 Verbrugge clamp 17.5cm
- 1 Hammer 200 G
- 1 Field type clamp 20cm
- 1 16 cm multitoothed tweezers
- 1 Stop end clamp 120 cm
- 2 Hexagon screwdriver tip with AO T8 anchor Ø2.5mm
- 1 Handle with torque AO green 0.8 NM
- 1 Osteotome mini-lamb Ø3 x 125mm
- 1 Osteotome mini-lamb Ø4 x 125mm
- 1 Osteotome mini-lamb Ø6 x 125mm
- 1 Osteotome mini-lamb Ø8 x 125mm
- 1 Watchmaker handle with blue AO input Ø2.5mm
- 2 Drill guide / sleeve Ø2.0 / 2.5mm
- 2 Drill bit with AO anchor Ø2.0 x 125mm
- 1 Lionter gauge 50mm
- 1 Ø 2.0/2.7mm Double Drill guide
- 1 Adson brown forceps
- 2 Mini-Hohmann spacer Ø 6mm
- 2 Kirschner Wire Ø1.1 X 100mm



internationalsales@fixier.com.mx
Mobile: (+52) 55 8106 3240

www.traufix.com

 traufix.mx

 traufix